

FASTKD2 Polyclonal Antibody Cat No: HR1AP3442

For research use only

Overview

| Product Name | FASTKD2 Polyclonal Antibody |
|-------------------------|---|
| Source | Rabbit |
| Applications | WB,IHC-p,ELISA |
| Species Reactivity | Human |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 81kDa |
| GeneID?Human? | FASTKD2 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | FASTKD2; KIAA0971; FAST kinase domain-containing protein 2 |
| Background | FAST kinase domains 2(FASTKD2) Homo sapiens This gene encodes a protein that is localized in the mitochondrial inner compartment and that may play a role in mitochondrial apoptosis. Nonsense mutations have been reported to result in cytochrome c oxidase deficiency. [provided by RefSeq, Oct 2008], |