



CYB561D2 Polyclonal Antibody

Cat No: HR1AP5699

For research use only

Overview

Product Name	CYB561D2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	24kDa
GeneID?Human?	CYB561D2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CYB561D2; 101F6; LUCA12.2; Cytochrome b561 domain-containing protein 2; Putative tumor suppressor protein 101F6
Background	cofactor: Binds 2 heme groups non-covalently.,similarity: Contains 1 cytochrome b561 domain.,