

## Tau (10E3) Mouse mAb

Cat No: HR1AM2432

For research use only

### Overview

|                       |   |
|-----------------------|---|
| Product Name          | Tau (10E3) Mouse mAb  |
| Source                | Mouse   |
| Applications          | WB, IHC   |
| Species Reactivity    | Human,Rat,Mouse   |
| Recommended Dilutions | WB 1:1,000-2,000 IHC 1:100-200  |
| Immunogen             |   |
| Species               | Mouse   |
| Storage               | PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.                       |
| Isotype               | IgG1  |
| Clonality             | Monoclonal  |
| Concentration         | 1mg/ml  |
| Observed band         | 50-80kDa  |
| GeneID?Human?         | 4137  |
| Human Swiss-Prot No.  |   |
| Cellular localization |   |
| Alternative Names     | Mixed lineage kinase domain like protein, FLJ34389  |
| Background            | Tau is a heterogeneous microtubule-associated protein that promotes and stabilizes microtubule assembly, especially in axons. |