

## Tau (10E3) Mouse mAb

Cat No: HR1AM2432

For research use only

## Overview

Product Name	Tau (10E3) Mouse mAb
Source	Mouse
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:100-200
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	50-80kDakDa
GeneID?Human?	4137
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Mixed lineage kinase domain like protein, FLJ34389
Background	Tau is a heterogeneous microtubule-associated protein that promotes and stabilizes microtubule assembly, especially in axons.