

IPP-2 (phospho Ser120/S121) Polyclonal Antibody

Cat No: HR1AP7807

For research use only

Overview

Product Name	IPP-2 (phospho Ser120/S121) Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	PPP1R2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PPP1R2; IPP2; Protein phosphatase inhibitor 2; IPP-2
Background	protein phosphatase 1 regulatory inhibitor subunit 2(PPP1R2) Homo sapiens Protein phosphatase-1 (PP1) is one of the main eukaryotic serine/threonine phosphatases. The protein encoded by this gene binds to the catalytic subunit of PP1, strongly inhibiting its activity. Ten related pseudogenes have been found throughout the human genome. Several splice variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2015].