



p130 Cas (phospho Tyr249) Polyclonal Antibody

Cat No: HR1AP9173

For research use only

Overview

Product Name	p130 Cas (phospho Tyr249) Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	130kDa
GeneID?Human?	BCAR1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	BCAR1; CAS; CASS1; CRKAS; Breast cancer anti-estrogen resistance protein 1; CRK-associated substrate; Cas scaffolding protein family member 1; p130cas
Background	BCAR1, Cas family scaffolding protein(BCAR1) Homo sapiens BCAR1, or CAS, is an Src (MIM 190090) family kinase substrate involved in various cellular events, including migration, survival, transformation, and invasion (Sawada et al., 2006 [PubMed 17129785]).[supplied by OMIM, May 2009],