

MMP2 (1H1) Mouse mAb Cat No: HR1AM2427

For research use only

Overview

Product Name	MMP2 (1H1) Mouse mAb
Source	Mouse
Applications	IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	IHC 1:200-500
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	64,72kDakDa
GeneID?Human?	4313
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	matrix metalloproteinase 2, MONA, gelatinase A
Background	The matrix metalloproteinases (MMPs) are a family of proteases that target many extracellular proteins including other proteases, growth factors, cell surface receptors, and adhesion molecules.