

# LRG1 Polyclonal Antibody

Cat No: HR1AP3836

For research use only

## Overview

Product Name	LRG1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	58kDa
GeneID?Human?	LRG1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	LRG1; LRG; Leucine-rich alpha-2-glycoprotein; LRG
Background	leucine rich alpha-2-glycoprotein 1(LRG1) Homo sapiens The leucine-rich repeat (LRR) family of proteins, including LRG1, have been shown to be involved in protein-protein interaction, signal transduction, and cell adhesion and development. LRG1 is expressed during granulocyte differentiation (O&apos;Donnell et al., 2002 [PubMed 12223515]).[supplied by OMIM, Mar 2008],