



## BTF Polyclonal Antibody

Cat No: HR1AP2908

For research use only

### Overview

Product Name	BTF Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	106kDa
GenelD?Human?	BCLAF1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	BCLAF1; BTF; KIAA0164; Bcl-2-associated transcription factor 1; Btf
Background	BCL2 associated transcription factor 1(BCLAF1) Homo sapiens This gene encodes a transcriptional repressor that interacts with several members of the BCL2 family of proteins. Overexpression of this protein induces apoptosis, which can be suppressed by co-expression of BCL2 proteins. The protein localizes to dot-like structures throughout the nucleus, and redistributes to a zone near the nuclear envelope in cells undergoing apoptosis. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008].