

# MAGE-C2 Polyclonal Antibody

**Cat No: HR1AP3852**

For research use only

## Overview

|                       |                                                                                                                                                                                                                                                                                                                          |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | MAGE-C2 Polyclonal Antibody                                                                                                                                                                                                                                                                                              |
| Source                | Rabbit                                                                                                                                                                                                                                                                                                                   |
| Applications          | WB,ELISA                                                                                                                                                                                                                                                                                                                 |
| Species Reactivity    | Human                                                                                                                                                                                                                                                                                                                    |
| Recommended Dilutions |                                                                                                                                                                                                                                                                                                                          |
| Immunogen             |                                                                                                                                                                                                                                                                                                                          |
| Species               | Rabbit                                                                                                                                                                                                                                                                                                                   |
| Storage               | -20°C/1 year                                                                                                                                                                                                                                                                                                             |
| Isotype               |                                                                                                                                                                                                                                                                                                                          |
| Clonality             |                                                                                                                                                                                                                                                                                                                          |
| Concentration         | 1 mg/ml                                                                                                                                                                                                                                                                                                                  |
| Observed band         | 41kDa                                                                                                                                                                                                                                                                                                                    |
| GeneID?Human?         | MAGEC2                                                                                                                                                                                                                                                                                                                   |
| Human Swiss-Prot No.  |                                                                                                                                                                                                                                                                                                                          |
| Cellular localization |                                                                                                                                                                                                                                                                                                                          |
| Alternative Names     | MAGEC2; HCA587; MAGEE1; Melanoma-associated antigen C2; Cancer/testis antigen 10; CT10; Hepatocellular carcinoma-associated antigen 587; MAGE-C2 antigen; MAGE-E1 antigen                                                                                                                                                |
| Background            | MAGE family member C2(MAGEC2) Homo sapiens This gene is a member of the MAGEC gene family. It is not expressed in normal tissues, except for testis, and is expressed in tumors of various histological types. This gene and the other MAGEC genes are clustered on chromosome Xq26-q27. [provided by RefSeq, Oct 2009], |