



Syndecan-4 (phospho Ser179) Polyclonal Antibody

Cat No: HR1AP8247

For research use only

Overview

Product Name	Syndecan-4 (phospho Ser179) Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	SDC4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SDC4; Syndecan-4; SYND4; Amphiglycan; Ryudocan core protein
Background	syndecan 4(SDC4) Homo sapiens The protein encoded by this gene is a transmembrane (type I) heparan sulfate proteoglycan that functions as a receptor in intracellular signaling. The encoded protein is found as a homodimer and is a member of the syndecan proteoglycan family. This gene is found on chromosome 20, while a pseudogene has been found on chromosome 22. [provided by RefSeq, Jul 2008],