




## Latexin Polyclonal Antibody

Cat No: HR1AP4975

For research use only

### Overview

Product Name	Latexin Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	25kDa
GenelD?Human?	LXN
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	LXN; Latexin; Endogenous carboxypeptidase inhibitor; ECI; Protein MUM; Tissue carboxypeptidase inhibitor; TCI
Background	latexin(LXN) Homo sapiens This gene encodes the only known protein inhibitor of zinc-dependent metallocarboxypeptidases. [provided by RefSeq, Oct 2008].



# bioelsa

