

HSP20 (phospho Ser16) Polyclonal Antibody

Cat No: HR1AP5830

For research use only

Overview

Product Name	HSP20 (phospho Ser16) Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	HSPB6
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	HSPB6; Heat shock protein beta-6; HspB6; Heat shock 20 kDa-like protein p20
Background	heat shock protein family B (small) member 6(HSPB6) Homo sapiens This locus encodes a heat shock protein. The encoded protein likely plays a role in smooth muscle relaxation. [provided by RefSeq, Jan 2012],