



Layilin Polyclonal Antibody

Cat No: HR1AP5063

For research use only

Overview

Product Name	Layilin Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	43kDa
GeneID?Human?	LAYN
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	LAYN; Layilin
Background	domain:The C-terminal domain interacts with the N-terminal domain of RDX.,function:Receptor for hyaluronate.,similarity:Contains 1 C-type lectin domain.,subcellular location:Colocalizes with TLN1 at the membrane ruffles.,subunit:Interacts with NF2, RDX and TLN1.,