

Crk-L (phospho Tyr207) Polyclonal Antibody

Cat No: HR1AP2628

For research use only

Overview

Product Name	Crk-L (phospho Tyr207) Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse,Rat,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	39kDa
GeneID?Human?	CRKL
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CRKL; Crk-like protein
Background	CRK like proto-oncogene, adaptor protein(CRKL) Homo sapiens This gene encodes a protein kinase containing SH2 and SH3 (src homology) domains which has been shown to activate the RAS and JUN kinase signaling pathways and transform fibroblasts in a RAS-dependent fashion. It is a substrate of the BCR-ABL tyrosine kinase, plays a role in fibroblast transformation by BCR-ABL, and may be oncogenic.[provided by RefSeq, Jan 2009],