

## YYAP1 Polyclonal Antibody

Cat No: HR1AP13069

For research use only

### Overview

Product Name	YYAP1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	87kDa
GeneID?Human?	YY1AP1 HCCA2 YY1AP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	YY1 associated protein 1(YY1AP1) Homo sapiens The encoded gene product presumably interacts with YY1 protein; however, its exact function is not known. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008],