

Pdcd-4 (phospho Ser67) Polyclonal Antibody

Cat No: HR1AP6660

For research use only

Overview

Product Name	Pdcd-4 (phospho Ser67) Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	56kDa
GeneID?Human?	PDCD4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PDCD4; H731; Programmed cell death protein 4; Neoplastic transformation inhibitor protein; Nuclear antigen H731-like; Protein 197/15a
Background	programmed cell death 4 (neoplastic transformation inhibitor)(PDCD4) Homo sapiens This gene is a tumor suppressor and encodes a protein that binds to the eukaryotic translation initiation factor 4A1 and inhibits its function by preventing RNA binding. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2010],