

TM171 Polyclonal Antibody Cat No: HR1AP12756

For research use only

Overview

Product Name	TM171 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	Dioeisa
Clonality	
Concentration	1 mg/ml
Observed band	35kDa
GeneID?Human?	TMEM171
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	TMEM171 (Transmembrane Protein 171) is a Protein Coding gene.