

## TM171 Polyclonal Antibody

Cat No: HR1AP12756

For research use only

### Overview

|                       |                                                               |
|-----------------------|---------------------------------------------------------------|
| Product Name          | TM171 Polyclonal Antibody                                     |
| Source                | Rabbit                                                        |
| Applications          | WB,ELISA                                                      |
| Species Reactivity    | Human,Rat,Mouse                                               |
| Recommended Dilutions |                                                               |
| Immunogen             |                                                               |
| Species               | Rabbit                                                        |
| Storage               | -20°C/1 year                                                  |
| Isotype               |                                                               |
| Clonality             |                                                               |
| Concentration         | 1 mg/ml                                                       |
| Observed band         | 35kDa                                                         |
| GeneID?Human?         | TMEM171                                                       |
| Human Swiss-Prot No.  |                                                               |
| Cellular localization |                                                               |
| Alternative Names     |                                                               |
| Background            | TMEM171 (Transmembrane Protein 171) is a Protein Coding gene. |