

BID (phospho Ser78) Polyclonal Antibody

Cat No: HR1AP8242

For research use only

Overview

Product Name	BID (phospho Ser78) Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	BID
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	BID; BH3-interacting domain death agonist; p22 BID; BID
Background	BH3 interacting domain death agonist(BID) Homo sapiens This gene encodes a death agonist that heterodimerizes with either agonist BAX or antagonist BCL2. The encoded protein is a member of the BCL-2 family of cell death regulators. It is a mediator of mitochondrial damage induced by caspase-8 (CASP8); CASP8 cleaves this encoded protein, and the COOH-terminal part translocates to mitochondria where it triggers cytochrome c release. Multiple alternatively spliced transcript variants have been found, but the full-length nature of some variants has not been defined. [provided by RefSeq, Jul 2008],