



TACC3 Polyclonal Antibody

Cat No: HR1AP4661

For research use only

Overview

Product Name	TACC3 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	80kDa
GenID?Human?	TACC3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TACC3; ERIC1; Transforming acidic coiled-coil-containing protein 3; ERIC-1
Background	transforming acidic coiled-coil containing protein 3(TACC3) <i>Homo sapiens</i> This gene encodes a member of the transforming acidic coiled-coil protein family. The encoded protein is a motor spindle protein that may play a role in stabilization of the mitotic spindle. This protein may also play a role in growth a differentiation of certain cancer cells. [provided by RefSeq, Nov 2011],