

Epsilon Tubulin (2B3) Mouse mAb

Cat No: HR1AM2414

For research use only

Overview

Product Name	Epsilon Tubulin (2B3) Mouse mAb
Source	Mouse
Applications	WB
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	55kDa
GeneID?Human?	51175
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	FLJ22589,TUBE1,tubulin epsilon 1, ? tubulin
Background	N/A