



## FASTKD1 Polyclonal Antibody

Cat No: HR1AP3441

For research use only

### Overview

Product Name	FASTKD1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	97kDa
GenelID?Human?	FASTKD1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	FASTKD1; KIAA1800; FAST kinase domain-containing protein 1
Background	caution:It is uncertain whether Met-1 or Met-15 is the initiator.,similarity:Belongs to the FAST kinase family.,similarity:Contains 1 RAP domain.,

