



## CIDE-B Polyclonal Antibody

Cat No: HR1AP3094

For research use only

### Overview

Product Name	CIDE-B Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	24kDa
GenelID?Human?	CIDE B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CIDE B; Cell death activator CIDE-B; Cell death-inducing DFFA-like effector B
Background	function:Activates apoptosis.,similarity:Contains 1 CIDE-N domain.,subunit:Inhibited by DFFB. Interacts with DFFA and DFFB.,tissue specificity:Highly expressed in liver and small intestine and, at lower levels, in colon, kidney and spleen.,