

CAGE-1 Polyclonal Antibody

Cat No: HR1AP2939

For research use only

Overview

Product Name	CAGE-1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	90kDa
GeneID?Human?	CAGE1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CAGE1; CTAG3; Cancer-associated gene 1 protein; CAGE-1; Cancer/testis antigen 3; CT3
Background	tissue specificity:Testis-specific expression in normal tissues, but wide expression among cancer tissues and cell lines.,