

TOB1 (phospho Ser164) Polyclonal Antibody

Cat No: HR1AP5538

For research use only

Overview

Product Name	TOB1 (phospho Ser164) Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	40kDa
GeneID?Human?	TOB1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TOB1; TOB; TROB1; Protein Tob1; Transducer of erbB-2 1
Background	transducer of ERBB2, 1(TOB1) Homo sapiens This gene encodes a member of the transducer of erbB-2 /B-cell translocation gene protein family. Members of this family are anti-proliferative factors that have the potential to regulate cell growth. The encoded protein may function as a tumor suppressor. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Aug 2011],