

## SLIK2 Polyclonal Antibody Cat No: HR1AP12382

For research use only

## **Overview**

Product Name	SLIK2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	92kDa
GeneID?Human?	SLITRK2 CXorf2 KIAA1854 SLITL1 UNQ9197/PRO34756
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	SLIT and NTRK like family member 2(SLITRK2) Homo sapiens This gene encodes an integral membrane protein that contains two N-terminal leucine-rich repeats domains and contains C-terminal regions similar to neurotrophin receptors. The encoded protein may play a role in modulating neurite activity. Alternatively spliced transcript variants encoding the same protein have been described.[provided by RefSeq, Feb 2010],