



# SLIK2 Polyclonal Antibody

Cat No: HR1AP12382

For research use only

## Overview

|                       |                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|-----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | SLIK2 Polyclonal Antibody                                                                                                                                                                                                                                                                                                                                                                                                        |
| Source                | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                           |
| Applications          | WB,ELISA                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Species Reactivity    | Human,Mouse                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Recommended Dilutions |                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Immunogen             |                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Species               | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                           |
| Storage               | -20°C/1 year                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Isotype               |                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Clonality             |                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Concentration         | 1 mg/ml                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Observed band         | 92kDa                                                                                                                                                                                                                                                                                                                                                                                                                            |
| GeneID?Human?         | SLITRK2 CXorf2 KIAA1854 SLITL1 UNQ9197/PRO34756                                                                                                                                                                                                                                                                                                                                                                                  |
| Human Swiss-Prot No.  |                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Cellular localization |                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Alternative Names     |                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Background            | SLIT and NTRK like family member 2(SLITRK2) Homo sapiens This gene encodes an integral membrane protein that contains two N-terminal leucine-rich repeats domains and contains C-terminal regions similar to neurotrophin receptors. The encoded protein may play a role in modulating neurite activity. Alternatively spliced transcript variants encoding the same protein have been described.[provided by RefSeq, Feb 2010], |