

DKKL1 Polyclonal Antibody

Cat No: HR1AP12304

For research use only

Overview

Product Name	DKKL1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	26kDa
GeneID?Human?	DKKL1 SGY1 UNQ735/PRO1429
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	dickkopf like acrosomal protein 1(DKKL1) Homo sapiens The dickkopf protein family interacts with the Wnt signaling pathway and its members are characterized by two conserved cysteine-rich domains. This gene encodes a secreted protein that has low sequence similarity to the dickkopf-3 protein. Multiple alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, Oct 2010],