

Muscle actin (3E9) Mouse mAb

Cat No: HR1AM2408

For research use only

Overview

Product Name	Muscle actin (3E9) Mouse mAb
Source	Mouse
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:100-200
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	42kDa
GeneID?Human?	58
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	N/A
Background	