



## BNIP-2 Polyclonal Antibody

Cat No: HR1AP2896

For research use only

### Overview

Product Name	BNIP-2 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	36kDa
GenID?Human?	BNIP2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	BNIP2; NIP2; BCL2/adenovirus E1B 19 kDa protein-interacting protein 2
Background	BCL2 interacting protein 2(BNIP2) Homo sapiens This gene is a member of the BCL2/adenovirus E1B 19 kd-interacting protein (BNIP) family. It interacts with the E1B 19 kDa protein, which protects cells from virally-induced cell death. The encoded protein also interacts with E1B 19 kDa-like sequences of BCL2, another apoptotic protector. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Mar 2016],