



## Bcl-2 (phospho Ser70) Polyclonal Antibody

Cat No: HR1AP2383

For research use only

### Overview

|                       |                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|-----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | Bcl-2 (phospho Ser70) Polyclonal Antibody                                                                                                                                                                                                                                                                                                                                                                                                         |
| Source                | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Applications          | WB,IHC-p,IP,ELISA                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Species Reactivity    | Human,Monkey                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Recommended Dilutions |                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Immunogen             |                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Species               | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Storage               | -20°C/1 year                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Isotype               |                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Clonality             |                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Concentration         | 1 mg/ml                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| Observed band         | 26kDa                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| GeneID?Human?         | BCL2                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Human Swiss-Prot No.  |                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Cellular localization |                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Alternative Names     | BCL2; Apoptosis regulator Bcl-2                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Background            | <p>BCL2, apoptosis regulator(BCL2) Homo sapiens This gene encodes an integral outer mitochondrial membrane protein that blocks the apoptotic death of some cells such as lymphocytes. Constitutive expression of BCL2, such as in the case of translocation of BCL2 to Ig heavy chain locus, is thought to be the cause of follicular lymphoma. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2016],</p> |