



PDIP2 Polyclonal Antibody

Cat No: HR1AP12104

For research use only

Overview

Product Name	PDIP2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	40kDa
GeneID?Human?	POLDIP2 PDIP38 POLD4 HSPC017
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	DNA polymerase delta interacting protein 2(POLDIP2) Homo sapiens This gene encodes a protein that interacts with the DNA polymerase delta p50 subunit, as well as with proliferating cell nuclear antigen. The encoded protein maybe play a role in the ability of the replication fork to bypass DNA lesions. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2014],