



## Fhit Polyclonal Antibody

Cat No: HR1AP6365

For research use only

### Overview

Product Name	Fhit Polyclonal Antibody
Source	Rabbit
Applications	IHC-p, WB, ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GenelD?Human?	FHIT
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	FHIT; Bis(5'-adenosyl)-triphosphatase; AP3A hydrolase; AP3Aase; Diadenosine 5'; 5'''-P1,P3-triphosphate hydrolase; Dinucleosidetriphosphatase; Fragile histidine triad protein
Background	fragile histidine triad(FHIT) Homo sapiens This gene, a member of the histidine triad gene family, encodes a diadenosine 5'&apost;5'&apost;&apost;5'&apost;5'&apost;&apost;-P1,P3-triphosphate hydrolase involved in purine metabolism. The gene encompasses the common fragile site FRA3B on chromosome 3, where carcinogen-induced damage can lead to translocations and aberrant transcripts of this gene. In fact, aberrant transcripts from this gene have been found in about half of all esophageal, stomach, and colon carcinomas. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Oct 2009],