

FOXF1 Polyclonal Antibody

Cat No: HR1AP11994

For research use only

Overview

Product Name	FOXF1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	41kDa
GeneID?Human?	FOXF1 FKHL5 FREAC1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	forkhead box F1(FOXF1) Homo sapiens This gene belongs to the forkhead family of transcription factors which is characterized by a distinct forkhead domain. The specific function of this gene has not yet been determined; however, it may play a role in the regulation of pulmonary genes as well as embryonic development. [provided by RefSeq, Jul 2008],