

BRMS-1 Polyclonal Antibody

Cat No: HR1AP4999

For research use only

Overview

Product Name	BRMS-1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	28kDa
GeneID?Human?	BRMS1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	BRMS1; Breast cancer metastasis-suppressor 1
Background	breast cancer metastasis suppressor 1(BRMS1) Homo sapiens This gene reduces the metastatic potential, but not the tumorigenicity, of human breast cancer and melanoma cell lines. The protein encoded by this gene localizes primarily to the nucleus and is a component of the mSin3a family of histone deacetylase complexes (HDAC). The protein contains two coiled-coil motifs and several imperfect leucine zipper motifs. Alternative splicing results in two transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008],