

ING4 Polyclonal Antibody

Cat No: HR1AP4878

For research use only

Overview

Product Name	ING4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	28kDa
GeneID?Human?	ING4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ING4; My036; Inhibitor of growth protein 4; p29ING4
Background	inhibitor of growth family member 4(ING4) Homo sapiens This gene encodes a tumor suppressor protein that contains a PHD-finger, which is a common motif in proteins involved in chromatin remodeling. This protein can bind TP53 and EP300/p300, a component of the histone acetyl transferase complex, suggesting its involvement in the TP53-dependent regulatory pathway. Multiple alternatively spliced transcript variants have been observed that encode distinct proteins. [provided by RefSeq, Jul 2008],