

KKCC1 Polyclonal Antibody

Cat No: HR1AP11919

For research use only

Overview

Product Name	KKCC1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	55kDa
GeneID?Human?	CAMKK1 CAMKKA
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	calcium/calmodulin dependent protein kinase kinase 1(CAMKK1) Homo sapiens The product of this gene belongs to the Serine/Threonine protein kinase family, and to the Ca(2+)/calmodulin-dependent protein kinase subfamily. This protein plays a role in the calcium/calmodulin-dependent (CaM) kinase cascade. Three transcript variants encoding two distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008].