

## LC3B (9H5) Mouse mAb

Cat No: HR1AM2402

For research use only

## Overview

Product Name	LC3B (9H5) Mouse mAb
Source	Mouse
Applications	IHC
Species Reactivity	Human
Recommended Dilutions	IHC 1:100-200
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	14,16kDakDa
GeneID?Human?	81631
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ATG8F antibody, MAP1 light chain 3 like protein 2 antibody, MAP1A/MAP1B LC3 B antibody
Background	Autophagy marker Light Chain 3 (LC3) was originally identified as a subunit of microtubule-associated proteins 1A and 1B (termed MAP1LC3).