

KT33A Polyclonal Antibody

Cat No: HR1AP11834

For research use only

Overview

Product Name	KT33A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	44kDa
GeneID?Human?	KRT33A HHA3-I HKA3A KRTHA3A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	keratin 33A(KRT33A) Homo sapiens This gene encodes a member of the keratin gene family. This gene is one of multiple type I hair keratin genes that are clustered in a region of chromosome 17q12-q21 and have the same direction of transcription. As a type I hair keratin, the encoded protein is an acidic protein which heterodimerizes with type II keratins to form hair and nails. There are two isoforms of this protein, encoded by two separate genes, keratin 33A and keratin 33B. [provided by RefSeq, May 2012],