

GSK3 β (2C6) Mouse mAb

Cat No: HR1AM2394

For research use only

Overview

Product Name	GSK3 β (2C6) Mouse mAb
Source	Mouse
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:100-200
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	46kDa
GeneID?Human?	2932
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Glycogen synthase kinase 3 beta antibody, GSK3B antibody, Serine/threonine protein kinase GSK3B antibody
Background	Glycogen synthase kinase 3 beta, also known as GSK3B, is an enzyme that in humans is encoded by the GSK3B gene.