

CD5L Polyclonal Antibody

Cat No: HR1AP11630

For research use only

Overview

Product Name	CD5L Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	38kDa
GeneID?Human?	CD5L API6 UNQ203/PRO229
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May play a role in the regulation of the immune system. Seems to play a role as an inhibitor of apoptosis.,similarity:Contains 3 SRCR domains.,subunit:Associated with IgM.,tissue specificity:Expressed in spleen, lymph node, thymus, bone marrow, and fetal liver, but not in non-lymphoid tissues.,