

UIMC1 Polyclonal Antibody

Cat No: HR1AP11538

For research use only

Overview

Product Name	UIMC1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	79kDa
GeneID?Human?	UIMC1 RAP80 RXRIP110
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	ubiquitin interaction motif containing 1(UIMC1) Homo sapiens This gene encodes a nuclear protein that interacts with Brca1 (breast cancer 1) in a complex to recognize and repair DNA lesions. This protein binds ubiquitinated lysine 63 of histone H2A and H2AX. This protein may also function as a repressor of transcription. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2015],