



## GNL3L Polyclonal Antibody

Cat No: HR1AP3553

For research use only

### Overview

|                       |   |
|-----------------------|---|
| Product Name          | GNL3L Polyclonal Antibody   |
| Source                | Rabbit  |
| Applications          | WB,IHC-p,ELISA  |
| Species Reactivity    | Human   |
| Recommended Dilutions |   |
| Immunogen             |   |
| Species               | Rabbit  |
| Storage               | -20°C/1 year  |
| Isotype               |   |
| Clonality             |   |
| Concentration         | 1 mg/ml   |
| Observed band         | 66kDa   |
| GenID?Human?          | GNL3L   |
| Human Swiss-Prot No.  |   |
| Cellular localization |   |
| Alternative Names     | GNL3L; Guanine nucleotide-binding protein-like 3-like protein   |
| Background            | G protein nucleolar 3 like(GNL3L) Homo sapiens The protein encoded by this gene appears to be a nucleolar GTPase that is essential for ribosomal pre-rRNA processing and cell proliferation. Two transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, May 2010], |