



GNL3L Polyclonal Antibody

Cat No: HR1AP3553

For research use only

Overview

Product Name	GNL3L Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	66kDa
GeneID?Human?	GNL3L
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GNL3L; Guanine nucleotide-binding protein-like 3-like protein
Background	G protein nucleolar 3 like(GNL3L) Homo sapiens The protein encoded by this gene appears to be a nucleolar GTPase that is essential for ribosomal pre-rRNA processing and cell proliferation. Two transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, May 2010],