



PNISR Polyclonal Antibody

Cat No: HR1AP11422

For research use only

Overview

Product Name	PNISR Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	88kDa
GeneID?Human?	PNISR C6orf111 SFRS18 SRRP130 HSPC261 HSPC306
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,sequence caution:Contaminating sequence. Potential poly-A sequence.,similarity:Belongs to the splicing factor SR family.,subunit:Interacts with PNN.,tissue specificity:Expressed in heart, skeletal muscle, thymus, spleen, kidney, liver, placenta and leukocytes.,