

Bax (6F11) Mouse mAb Cat No: HR1AM2379

For research use only

Overview

Product Name	Bax (6F11) Mouse mAb
Source	Mouse
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:100-200
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	20kDakDa
GeneID?Human?	581
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Apoptosis regulator bax, BAXA, BCL2 associated X protein
Background	Bax is a key component for cellular induced apoptosis through mitochondrial stress. Upon apoptotic stimulation, Bax forms oligomers and translocates from the cytosol to the mitochondrial membrane.