

SPB10 Polyclonal Antibody

Cat No: HR1AP11345

For research use only

Overview

Product Name	SPB10 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	43kDa
GeneID?Human?	SERPINB10 PI10
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	serpin family B member 10(SERPINB10) Homo sapiens This gene is a member of the serpin peptidase inhibitor, clade B family and is found in a cluster of other similar genes on chromosome 18. The protein encoded by this gene appears to help control the regulation of protease functions during hematopoiesis. Variations in this gene may increase the risk of prostate cancer. [provided by RefSeq, Dec 2015],