



## Exo1 Polyclonal Antibody

Cat No: HR1AP3420

For research use only

### Overview

Product Name	Exo1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	94kDa
GenelD?Human?	EXO1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	EXO1; EXO1; HEX1; Exonuclease 1; hExo1; Exonuclease I; hExo1
Background	exonuclease 1(EXO1) Homo sapiens This gene encodes a protein with 5' to 3' exonuclease activity as well as an RNase H activity. It is similar to the <i>Saccharomyces cerevisiae</i> protein Exo1 which interacts with Msh2 and which is involved in mismatch repair and recombination. Alternative splicing of this gene results in three transcript variants encoding two different isoforms. [provided by RefSeq, Jul 2008].