



Exo1 Polyclonal Antibody

Cat No: HR1AP3420

For research use only

Overview

Product Name	Exo1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	94kDa
GeneID?Human?	EXO1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	EXO1; EXOI; HEX1; Exonuclease 1; hExo1; Exonuclease I; hExoI
Background	exonuclease 1(EXO1) Homo sapiens This gene encodes a protein with 5' to 3' exonuclease activity as well as an RNase H activity. It is similar to the Saccharomyces cerevisiae protein Exo1 which interacts with Msh2 and which is involved in mismatch repair and recombination. Alternative splicing of this gene results in three transcript variants encoding two different isoforms. [provided by RefSeq, Jul 2008],